

Product datasheet for SC403229

MIR218a-2 Rat MicroRNA Expression Plasmid (MI0000957)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GACCAGUUGCCGCGGGCUUUUUUGUGCUUGAUCUAACCAUGUGGUGGAACGAUGGAAACGGAACAUGGUUCUGUCAAGCACCGCGAAAGCAUCGCUCUCUCU
Mature Sequence:	UUGUGCUUGAUCUAACCAUGU
Full insert:	>MI0000957 TATATTATGGGACAGATTTTACAGGCTCAATCCTGGTTCAGTATGTAAGTGGTGCCTCAGCTGCCAGGAT GTTTCATGGATATGAGTATCTGAAGGTTTCCACCTTGTCCCTCTCCTGCTGAGTCTCCCTGGGGCAGG GCACTGGGAGCGGGGCCCTGTCTTTAGCAGGATCCAGTGGACCGCACTACTGGTACTAATCACGTTCA GTGGGGCCTGCTCTGGCTTCTGCTTCCACGCTGCTTCTCTGAGCGCTCCTGTCTCTCTGACGC TGCTTCTGACCTTGACTCTGACCAGTTGCCGCGGGCTTTCCTTTGTGCTTGATCTAACCATGTGGTGG AACGATGGAACGGAACATGGTTCTGTCAAGCACCGCGAAAGCATCGCTCTCTCCTGCAGCATGGCCCA CCGCCCGCCGCCACCACCACCACCCTGGACTCTTCTCTGCTCTGGATCACCGCGCTGCCAGACCCA CTTCTCCCAGTCTTAACTCCCAGAGGGAAGGGTGGGTGGGAATCCTGAAGGAGCCAGGTACATGGTGAC AGCAGTTTAGAGAAACCCACGGGGCCATGCTAAGAGTATCCACATAACTTAATGAAATGTAGTCACCAT GCTCAATTCTGGTCTTGTGGGACCCAGGAGTGACAGCTAGGTTCCGACTTCCATTTGAAGAAGACA AGTGTACCTG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000957
Synonyms:	Mir218a-2 rno-mir-218-2 rno-mir-218a-2



[View online »](#)

Locus ID: 100314057