

Product datasheet for SC403230

MIR218b Rat MicroRNA Expression Plasmid (MI0015428)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGUAGAAAGCUGCGUGACGUUCCAUGGUGCUUGACGGAACCAUGUUUUACUCAUACCUCGCAAGCACAUUGGUUAGAUAAGCACAAACAGAAAAUCUCGCUACGUU
Mature Sequence:	CAUGGUUAGAUAAGCACAA
Full insert:	>MI0015428 TAGTCATCTATTTACAGATAATCACACGATGGGTATGAATTGTTAAAAAAGGTTTATTTTCTCTG AGTCAAAAATTGAAAACATATCTGGTAATAAACATACTCCTGCATATTATTGTAAGTGGGGCTTTCAAG GTTTTAAATTTTTTCTACAGACTTTTGTATTCTTCCCTCTCACATATTCTTCAGAAACGCCCCAAA GCCCTCTTTAATTTTTTTTTCATTTTGAGTCTATAAAAAAAAAAAAAAGACAAGTAAGTCCTGCCTTGA GCCTCAGCAGCTTGTCATGCTGTAGAAAGCTGCGTGACGTTCCATGGTGCTTGACGGAACCATGTTTTAC TCATACCTCGCAAGCACATGGTTAGATCAAGCACAAACAGAAAATCTCGCTACGTTATCACTACACCAACC ACTAGGCCAGAGGGAAGACCGAGAGGAAAGACAAAGGAGGAGAAAAGCTACGGACGGTCCGGAAGCATA GTACGGTTACGTGACATGCATTCTTCCCGCAGATTGTGTATGATCATAGCGGAAGCATAACTACCAAAGG TGCCTGGAAACACAAAGAAGGGGCTATGTGGAAAAGGCGCCTCCCATGTTTCTCAGAACTGCAGCGTGGG CAACCACAGCTTTCGGGAAACATCCAGGAAGAAGTTGAAGGTGAGCTACAAGGACAACAGAGAAGCTCAA GAGATCTGTGGTTT
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
Disclaimer:	flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015428
Synonyms:	Mir218b rno-mir-218b rno-mir-3565


[View online »](#)

Locus ID: 100526634