

Product datasheet for **SC403240**

MIR23b Rat MicroRNA Expression Plasmid (MI0000853)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Rat

Vector: pCMV-MIR (PCMVIR)

Precursor: CUCACCUGCUCUGGCUGCUUGGGUUCUGGCAUGCUGAUUUGUGACUUGAGAUUAAAAUCACAUUGCCAGGGAUACCACGCAACCAUGACCUUGGC

Mature Sequence: GGGUUCUGGCAUGCUGAUUU

Full insert: >MI0000853 ACCATGCCTACTGAGAAGAGTTTCCTATCCAAGGGGCCATAGACATCTCAGGCCATTAATAATCTGCAG
TTTGAGAACAGGGTGTGTCCCGAGGTTAAAGGTCACATACTCAGAGCGAGGATTGGTCATGTAGCTGAT
AAGTGTGCCTTTGTCTCCTGGTCCCTATGCTCCTCTATACAGTGCAGCTAGCATGCTGGATGACAGCA
GGCTACATACACTGCTAGCGTGGTGGGCTCAGAGCGGGTGGAGGACACTCAGCGCATGTGGATGGGAG
TGGTTGGGCATGGCGCTGCCCTCACCTGCTCTGGCTGCTTGGGTTCTGGCATGCTGATTTGTGACTTGA
GATTAATAATCACATTGCCAGGGATTACCACGCAACCATGACCTTGGCTGCTCCTGCAGAACCGTGGCTAC
TTTCTCTCTGCAGTTTGGAGAACAGAGGCACCTTGTGGCTCTTTGGGAAAAAAAAAAAAAAAAAAGAAGC
CACCAGCCTTCGAAGATGCTCACCAGCCCTTAATTATGGCCGGCGATGACCTCTTAACAAGGTGCAGA
GCTTAGCTGATTGGTGAACAGTATTGGTTCCGCTTTGTTACAGTGGCTAAGTTCTGCACCTGAAGAG
AAGGTGAGATGGGACAGTTAATTGCAGCTGCCAGGCAGAGGCCACTGCCTGGTGGTCACTCCAGA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000853](#)

Synonyms: rno-mir-23b

Locus ID: 100314002



[View online »](#)