

Product datasheet for **SC403253**

MIR293 Rat MicroRNA Expression Plasmid (MI0012609)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVIR)
Precursor:	CUGUGAUACUCAAAACUGUGACACUUUGUUUCUUUUGGAAGAAGUGCCGCAAAGUUUGCAGUGUUACCGAU
Mature Sequence:	ACUCAACUGUGACACUUU
Full insert:	>MI0012609 AAAGAATTGGACAACCCAGGGTGGATATGGTGATTTATCTACCCATGCACACTTCAATGGGGATGTTCTCA GTAGGGTAGCCTCCCATGAGGGCTGGTGCAGTGTTCCTGGAGCGGAGGACCGAGGGCCTTAACCGCT GAGCTAAATCCCAACCCACAGTGAATTTTCGTTGTCCTTGCTTTTCAGGTTTGTATAGGATTCCAG TGAGTTGAACTGTATTCTACTTAGCCTGGGCTACTCTACAACCTGCATGCCTGGCTATAAGCTAGAAAAC TTATTTGCAAGTGTTCACCTGTGATACTCAAAGTGTGTGACACTTTGTTTCTTTTGGAAAGTGGCCG CAAAGTTTGCAGTGTACCGATGGAGAAAAGTGCAGTTAGGGGGCTGTGTGGCTACTCACCTGGTGGCAC TCAAATGTAGGAGCGCTCTTCTGATCTGGTTAAGTGGAAAGTGTGCATCGTTGGGTGTCATCTGTGG GTTATGATTTCCACAAGTGCAGTCCGGGAAAGTTGAGTTGGTAGGGGAAAAGAGTACTGAGCCTTCTGAG ACTGGAAGGATTTGACAGGAGTACTGTATTTAACTGCAAGTGTTCATGAGATTGAGAGATTGAATTGC AGGCATGAATGGAGGATCTTGGAAAAACCACCATAGCCATGC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0012609
Synonyms:	rno-mir-293
Locus ID:	100314136



[View online »](#)