

Product datasheet for **SC403277**

MIR30a Rat MicroRNA Expression Plasmid (MI0000870)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVIR)
Precursor:	GCAACUGUAAACAUCCUCGACUGGAAGCUGUGAAGCCACAAAUGGGCUUUCAGUCGGAUGUUUGCAGCUGC
Mature Sequence:	UGUAAACAUCCUCGACUGGAAG
Full insert:	>MI0000870 CATGTTGTAGTCCTAGTAAGTCACCTCACCCTTGCATGATTGTTTCTTGTGGCCAACAACAGTAAGCCA TTTACATGTGTTGTTTTTCATAATATGTTTGAGATAGATAAATAATAATAAATGTTCTGAAACATATAAA CTATGTCTTATATAAGCTCTAGTTATATTCTACTAGGGCATATCTGAACGAGGCTTTACAGTTTACAGAA TGTTGCCATTACATTTTAGAAACACCTACAGGGGACTTCTCAGGTCAGTAAGACAGAATGCTGAAGAAGA CATATCAATGTTGACAGTGAGCAACTGTAACATCCTCGACTGGAAGCTGTGAAGCCACAAAATGGGCTTT CAGTCGGATGTTTGCAGCTGCCTACTGTCTCCTCAACTTCAAGGAGCAATTATCTTGTGGCTGAAACCAAT ACCTTGCTTTCTCTTTGATTTATTTTTACAAAGCATGGATGGAAAGAATAAACAACCTTCAAACCCACAGA AATGTTCTCCTGTGCACAGAATATGCAATATATGTACCTTAAAAATATAACAGTATTCTTAAGATAAGC AAAGAAATGTGTTATAGCATTGTTAGTTTTTAAAAAATAGGACAGTGCTTCAAAAAGCTATTTTTAAAA AGTAGGCTCTTAATCCACTATTTTTAAATCGAAAATAGAT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000870
Synonyms:	rno-mir-30a
Locus ID:	100314011



[View online »](#)