

Product datasheet for **SC403283**

MIR3120 Rat MicroRNA Expression Plasmid (MI0015476)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: GGUUGUCAUUCAGGCCUGGGUUGUCAUGUGACUGCCUGUCUGGCCUGCUGUACAGGUGAGCGGAUGUUCUGCACAGCAAGUGUAGACAGGCAGACACAUGACAAC

Mature Sequence: CACAGCAAGUGUAGACAGGCA

Full insert: >MI0015476 TTTTCATGTCAGTGTCTAGGAGGCAGCAAAGGTGTGTACAATGAATGGTCAGTGTGGTAGAGGAAAAGA
GGACCACGTCAGACACAGGCACATTTTCAGAACTGGCTTTATTAGGTTTCATAGGCACCACTCCCTTTA
CTTGGGGGTAGGTTTTAATATACGAATCAAATAGAAGGATTACATACGAGTTTATTTAATAAGCCGA
CAGAAAAGCCCAGGCCCCCAACCCCTCACTGTGGTTGTAGCTCTTGGTGTAGATGCTATGGTGTGAGG
GCAGCTTTCTTCAATGGCTGGTTGTCATTCAGGCTGGTTGTCATGTGACTGCCTGTCTGTGCCTGCTG
TACAGGTGAGCGGATGTTCTGCACAGCAAGTGTAGACAGGCAGACACATGACAACCTGTCCATCCAGGA
CCATGGCTCCTTCGCGTTCCTTTTTCTTGGACAGATACAATCAGCCTGTTTTCCAGATTTGGTGATTG
AAGGAGAGGGGATAAAGGAAAGGAGAAAAACACACAATAGAATTAACATTAAACATTATGGAATAAA
AAGAAGTTAGACGTTTACTTTATTACAGGAAGTACACTCATAACGTCATCAGAATGTTGCTAACAAGGA
TGTTGTCAAAAACAGAAAATTCACACGATTTAAGGAAAAAAAACAGCATAATTTTTTTAAAAATCC
AACTCCTATCGTTACTG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should
affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015476](#)

Synonyms: rno-mir-3120



[View online »](#)

Locus ID: 100526611