

Product datasheet for **SC403291**

MIR325 Rat MicroRNA Expression Plasmid (MI0000596)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: AUAUAGUGCUUGGUUCCUAGUAGGUGCUCAGUAAGUGUUUGUGACUAUUUCGUUUUAUUGAGCACCUCCUAUCAAUCAAGCACUGUGCUAGGCUCUGG

Mature Sequence: CCUAGUAGGUGCUCAGUAAGUGU

Full insert: >MI0000596 ATAGGTTTGCAAGTGTTAAGATAGGGCTGTTCCAGAAAGCCTTCTGTCAGCTATCAGATCCCATATCAAAC
ATACACATGACCATTTATATGAGCAAATCTCAAATTAGTCTTTTGGATTCTGCTGGGTTTTATCAGGCT
TATTCACAAAATTCAAGTGCCTGTTGTAAGGAGGGTCCAGAACAGAGCAGTAGTTGTGGGCATTCACAA
AGAGAGAAGGTTTTCTGCGAAAACCAATGGTTCCTGAATGAGGGTTTTGTGATGGATTAGGGGCCAAGTGT
TGTGGATTCTAACAGAGCCGATATAGTGCTTGGTTCCTAGTAGGTGCTCAGTAAGTGTGTTGTGACATAAT
TCGTTTATTGAGCACCTCCTATCAATCAAGCACTGTGCTAGGCTCTGGGGCCACATAGCTAAGTTAGACA
AGATGTTCCCTTTTGGGTCTTATTTCTAGTAGTGGCTATAATAGAATGATTAGTTATAACAACAGTGC
TGAGAGGGATAAGAGGATATGCCAGGAAACAGAAGCAAGTTGCACTGTGGGAAGAAGCTGCGAATATAC
ACCATTTATTCCAAAACCTTAGGAAGGAAAGGAGAGAAGGTGACAGCAAATAGAGAGGTAAGGGTGAGAA
GTTTTACCTTGAGAATGAGAGGCTTGAGCATAATGTTCCATTAGAGAGGGCAGGATGAAGATTCAAG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000596](#)

Synonyms: rno-mir-325

Locus ID: 100313970



[View online »](#)