

Product datasheet for **SC403306**

MIR344a-1 Rat MicroRNA Expression Plasmid (MI0000629)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: CUGCAGCCAGAGUUUUUACCAGUCAGGCUCCUGGCUAGAUUCCAGGUACCAACUGGUACCUGAUCUAGCCAAAGCCUGACCGUAAGCUGCAAAGAAA

Mature Sequence: UCAGGCUCCUGGCUAGAUUCCAGG

Full insert: >MI0000629 AGTACAACGGCAGTTTTTGTGTATTGTTAATTAGAGCCAGTGCTCAGGGTTAAATTTTCATTTTTTTGCT
GCCTTTTTCTTTTCATTTTCCCTGCTTCTGCTACAAGACAAAGGCAGGGGGTGGAGGGTGGGTAATGAA
TCAGCAGCTCTTGTGGAGATCTCCAAGGTTACCGCAGGCTTGCAGGCTTTTGCCTTCAAGCCTAGGTTT
GAGATCTCTCTCCTGATTGGCTCTGAGTCCATCCTTATTGAATGCACGCACTCCTCCAATCAGAGGCA
AGCAGAGGTAATAAAAATCTCTGCAGCCAGAGTTTTTACCAGTCAGGCTCCTGGCTAGATTCCAGGTACC
AACTGGTACCTGATCTAGCCAAAGCCTGACCGTAAGCTGCAAAAAGAAAAGGAGGACTGCTCCTGAACATT
GCAGGAGACAGAGCCAATGAAACAGCATGGGCAGATTCTTCTGGGATCCCCAGTCTCCTGCTTTAGTCAC
AGGCTCCTGCTTCTGTGTAGGAAATGTGAATCTGCTTCTGAAGCCATTGCTTCTGCTTTGATATCCAGC
ATCCTGTCTCAGTCTGTAGCCCTGCATTAATAACAAGGAATGTGATAGAGTTTGAATCCATTTCTGA
GAGAGATTATGGAAGTGCTGTTCCCTCACATTTCTTGGGGACAATACTGTGAAAACTGCTATGGCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000629](#)

Synonyms: Mir344a-1 rno-mir-344-1 rno-mir-344a-1

Locus ID: 100313981



[View online »](#)