

## Product datasheet for **SC403318**

### MIR34c Rat MicroRNA Expression Plasmid (MI0000876)

#### Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVIR)
Precursor:	AGUCUAGUUACUAGGCAGUGUAGUUAGCUGAUUGC UAAUAGUACCAAUCUAACCACACAGCCAGGUAAAAAGACU
Mature Sequence:	AGGCAGUGUAGUUAGCUGAUUGC
Full insert:	>MI0000876 GGCTTCCTTTAGGGACCGTTCGGCTGGCTACCTGTGTGTACGAGCCAGGGGAAAAAAAAATAACGTTTCG GTGATGTGTAAGAGGATCGGTATACACCGCCTTAACATTA AAAAAAAAAAGGATGCCAGGAAGAGCCCC CCAGAACAGTTCCTGCTGCGTGTGATTCACTGCCTGTTTTCCACTTCTGAGAAAGCTGTCAACCAAGG CAGCGACTAGAGTCAACCACTGCTTGGCTGAGGATTTGCTGGTCAACCACTGAGAGCACAACCAAG ACTGCAGGCTTTTTTCTATGAGTCTAGTTACTAGGCAGTGTAGTTAGCTGATTGCTAATAGTACCAATCA CTAACACACAGCCAGGTA AAAAGACTGAAATTCCTCCGACAGAGCCTCCTGTGCATCATTTTGGTGAGG ACAGCCACAGGGAAATAGCCTTACCATCTACTACCTATTAGTAGTTTAGGAGAAAAAACACGCTAGCCC TTAGCTGCTAAGTATACAATGTATCAATATGTGTGGGTACATTTTTGGTGATTTTTTAGATATATTGGA AAATAGAAATAATTACCCCTTAGCTGATTTGGTCATTATTGACTTGTGATTATAATCTCTGGACTTT ACCTGACACTATATTTTATTAATAACTCTAAAAATAAAGCGCCA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	<a href="#">MI0000876</a>
Synonyms:	rno-mir-34c
Locus ID:	100314014



[View online »](#)