

Product datasheet for **SC403334**

MIR3556b Rat MicroRNA Expression Plasmid (MI0015414)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: GCUGGUGCUUUUCCCCUACAUCAUAACCGAUUUCAAUUGGUGCUAGACAAAAACAGACUCUGAACACCAGGAGAAAUCGGUCAGCCUGUGUAAGAGAUGGACCA

Mature Sequence: ACCGAUUUCAAUUGGUGCUA

Full insert: >MI0015414 AATAGGTAAAACACAGACATGATAGATTGAATGTAGCAAATTATTTATTAAGCCAGCACTTAACAGACG
ATTCAACAGCAACTAGTTTATCACGGTACTCCATGTCAGAAATTGTAAGTACTGAGCATATTTACAACAATA
AATTACAACAATTAACAAAAAACCAGTAAGTGATCTCTAACTTCTGTAGCCACATACATATGCATTA
AAGAAGTTAGGAACTGTGATCTTCTTTCTGCTCTGCCTGTGTAGGAAATGACAGGCTGATGCTGACC
AGACACCACAAGGCTCTGAGGCTGGTCTTTTCCCCCTACATCATAACCGATTTCAAATGGTGCTAGACA
AAAACAGACTCTGAACACCAGGAGAAATCGGTACGCCTGTGTAAAGAGATGGACCAGTGGTGTCTTCAGG
TCTGAATCTCAAGGCAGGTGCCGGTCTGATTTCAAGCTGTGCTGCAACCCAGGACTGAGTGTCCACGGATC
CTCCAGGTACAGGCAGGCACTGATGTGAATGAAGCAGTCCACATACAGAGCCCTTACCCTCTTC
CTTGCTACAATTGACACCTGGAACATTACCAGACAATGAACAGTGCAGGCCAGAGGCCACCTCTCCACC
CATGGAATGCTGGGTAAGTGGATCACTCCGTACCCCTGATCACACCCTTACATTTATTCTGCACAAAA
CTTTCTC

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#).

Disclaimer: The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015414](#)

Synonyms: rno-mir-3556b



[View online »](#)

Locus ID: 100526572