

Product datasheet for **SC403347**

MIR3572 Rat MicroRNA Expression Plasmid (MI0015442)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGCCACCAGUAAGUGGAGGCCAGUUUGGGGGACAGGGCAAGGUGGGCAGUAGCUGACAGCCUGUUUACACUUGCCUUUUUCCCCAGGCUAUAAGUUCGGUUGUG
Mature Sequence:	UUACACUUGCCUUUUUCCCCAG
Full insert:	>MI0015442 CCCCTCTGGCCCCCACCTGCCACTTCTGCCTTGAGCTCCACCCTGGCAACAATGGGGCCTCTACCCCA GCAGCACCTCCAAGTGCCTGGGCCCTAAGGCCAGTCCAGCTCCAGCCACAACCTCGGGTACTCCTGCCC CCTATGCCAGGCTGTGGCCCCACCAATGCCAGCGGGCCAGCAATGCCAGCCCCGGCCCCCAGTGC CCAGCCAAGCGGGGAAGTGGTGGTGGCAGCGGAGGAGCAGCAGCAATAGTAACAGTGGCACAGGCGGA GGGGCTGGCAAGCAGAACGGTCCACCAGTAAGTGGAGGCCAGTTTGGGGACAGGGCAAGGTGGGCAGT AGCTGACAGCCTGTTTACACTTGCCCTTTTTTCCCCAGGCTATAGTTCGGTTGTGGCAGACAGCCCTGCA GAGGTGGCCCTGAGTAGCAGTGGGGGAGCAGTGCCAGCAGCCAAGCCCTGGGGCCACATCGGGCCCTC ATAACCAGCTCCCAGCACCTGTGAGTTGTGGGGTAACAAGTGTGTGGGCACAGGCTCTCATGCCACTC TTGTTTTGGGTCCAAGTTCGGCCCTTCTTGTGTAGTATCCAGTGTATCCTTTACCTCCACAAGGGTTC TTGTCCTGCTTCTTAGGGACTAGGACATGTGTACATTAGTCCCAGCTTGTCTTTTTGTTTTTTACT GTGGGTTT
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
Disclaimer:	flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015442
Synonyms:	rno-mir-3572



[View online »](#)

Locus ID: 100526587