

Product datasheet for **SC403355**

MIR3580 Rat MicroRNA Expression Plasmid (MI0015453)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMVMIR)
Precursor: GUGUGCUAUGCCUGCUUGUACUUAUUCUGGUGACAUGAGGUAUCUGUGACUAGGGUAGUAUGAGUAGAGUAAAGUUCACC
Mature Sequence: UGCUUGUACUUAUUCUGGUGAC

Full insert: >MI0015453 GAAAGAGAAGAAGTACTTGGGAGAGCTTTTCTCAAAATGAAAAACAAGACTAAAGAAACATCAATGTAA
TTCATAAATGCAAATCTGCCTAAGGCTTGTGTCTTTGCCAGTCAGGCTTTAAAAGAGAATGCAAGTAGA
TGTGATTTTAGTAAAAGTATAATGAAACAATTAGATGACTATTTTCATCAAGAAATACATCAAAAAGGG
GGTATGATTTTTTATTGGGATAATTCGTTTTATTTTGAATATTGAAAATGATGGAGGGAGTGTGATGGG
ACTTCCTCATTAAATGTAATAGTGTGCTATGCCCTGCTTGTACTTATTCTGGTGACATGAGGTATCTGTGA
CTAGGGTAGTATGAGTAGAGTAAAGTTCACCAAATCTTCATTTGATCAGTCTGGATGATGAGTATGGTG
GCTTCTCCTTAATTTAGATGTTCTATGGTTACTTCTGTAACCTAGGCTTCCTTTGGATGATTTTCAGTTGT
TCATTTGCCATCTCACAATAGCTAGAAGTCAGAGAATTGAACACAGTTATTTAAAATGTTTTATAATAGA
TAAATATTTAATGTGATCTCATTATGCTTTGACAATTATAGAAGCAAACAGAACTTGACTATAAAATTC
AGGTGGAAAAAATGCAAGATGGAAGACATCATACCAGAGTTGTATGATTGT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0015453](#)
Synonyms: rno-mir-3580
Locus ID: 100526639



[View online »](#)