

Product datasheet for **SC403359**

MIR3586 Rat MicroRNA Expression Plasmid (MI0015461)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CGGGAGGUGAGUAUGACGAUCCUGUCCUCAAGGAGCCUCAGUCUAGUAGGGAGGAGACAUACUAGACUGAGCUCCUCGAGGGCAGGAAAGCGCUCAUCAUC
Mature Sequence:	UCCUGUCCUCAAGGAG
Full insert:	>MI0015461 TTGAGGGGACAAAGTTGGAATCATGTTTGATTTCTTGGTGTGGCAAAGTGCTACTTGCTTCTTAGTAA AGAATTGGATGTAATAAAACGAAGCAAACCTCAATAACAGAACTCATTTAATTATTGTGCATGTTGTTT ATGGAAAACCAACAAGGAATTGTGTTATAAAGTACACAGTAAAATGTCTTATGTTGGCCATACTACAAA TGTTCCAAAGGTGAAGTCCAACCTCCACTGCAGTCACTGAGTACAAGATAGCCATGGCTACAGTAGAAGCA CATGTGCACTGTGCTCCTGGCGGGAGGTGAGTATGACGATCCCTGTCCTCAAGGAGCCTCAGTCTAGTAG GGAGGAGACATACTAGACTGTGAGCTCCTCGAGGGCAGGAAAGCGCTCCATCATCTGTGCTCCATGCCCA GGACCTGAGTAGCCAGGTGCAGCTGGCTGGAGCATGGCTAACCCAAGAGAGCTACAGGGACTTGAAGAA CAGGCTGCAAAGTCACAGAAGGGCCAGCCAGGACACAGAGCCAAGAAAAATGGGATTCTGCCAGTGGAAAG ATGCATTTTATGTGTTTAGAGTCTGGCCTTGTTTAAGAGGTCAAGATATTCTTGGGCCCTAAATGTATCC TTTCAGCTTTTGCCTGGACACAGAAAGAAGCCAGTGGAAAGGTGCTTGAGCAGCATCAATACTTTCTAGA AAATCATACTGAACCAAGTGA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015461
Synonyms:	rno-mir-3586



[View online »](#)

Locus ID: 100526641