

Product datasheet for **SC403360**

MIR3587 Rat MicroRNA Expression Plasmid (MI0015463)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMV MIR)
Precursor: UUGUGUGGUCGGCAGCGAAGUUGUUCUCCACCAAGAACACUGAUUUCAAUGGUGCUAGACAA
Mature Sequence: AACACUGAUUUCAAUGGUGC
Full insert: >MI0015463 TGAACACCAAAAAGAAATCAGTCATCCTCTAAGGGGTACAAGAGGTCATGTGCATTGGGAAGGCCCGGGC
CTTCTGTCTGTTGTACATGTAGGGTCTTCATCCAGGCTTGATGGAGTCTGCTTCCAGAGTAGGCTTACTA
TGTGAACAAGTATGAAATTCTGACAACGGATGACTAATTAGCTATCCACAACACCGAGGTAAGCAGGGCG
GGCAGGGCAGGCCCGCTGCTACTTGAGCAGGTTCTGGTGGCCTGACTGTCTTAGATACTGAATACCTCA
GGGCTGTCCACTCCTTCTTTGTGGTGGCAGCGAAGTTGTTCTCCACCAAGAACACTGATTTCAA
TGGTGCTAGACAATCACTATTTAAATCTAAACCACCATATGAAACCAGCTTCTGAAGAAGCTTTGTCTG
TTCTCGTCTCTTCCAACGGTCCCCTACAGAGGGGAAAGGTGGCGCTCCTATGCTGGCCTCATCCAGCAG
GTGGGGGGGGGGTGACATGTCAGTCAGTCACTTGTGCATTCTCCCTCCACATGATAACTGACTGGGCCCTC
ACCTGAGAGCAGCCGTGGGACATTGAGCCGTGACGCTTGACATTCACTGTGGCAGAGCTGAGCACACTA
CAAGCATGGCTTCTAGATGAGTCTGTTAGACTC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015463](#)

Synonyms: rno-mir-3587

Locus ID: 100526602



[View online »](#)