

Product datasheet for **SC403370**

MIR3596b Rat MicroRNA Expression Plasmid (MI0015455)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: GGCCCUAAGAAAGGCAGCAGGUCGUUAGUUAAACUCCUUGUGGGCAAAAUCUCUGCCCUAAAACUAUGCAACCUACUACCUCUCCUAGGAGCCC

Mature Sequence: AACUAUGCAACCUACUACCU

Full insert: >MI0015455 GGAAATAGCTGATGCTCCCATGTCAAAGACAAGGAAAATTAAGTAAAAATTAATATGAAAATTCTGCT
TTCCCAAACCTACCCAGTAACCCCAAATGTGAGACACAGGAGCAACTAGGGTTAACTTTGGAACCTCAG
ATTCTTACTAGATGCCTGTCGCTGACGCAGTGCAGTGTGCTGTGTACTAACCAGTGCACACTGTCTT
CAAAAGGTCTGCAATGTTCCATTCTTACGAGAACCACCTGCACGAAGAACTCTTGCTTACACGCATGGA
ACTGTATGCTGGTATGATAAGGCCCTAAGAAAGGCAGCAGGTCGTATAGTTAACTCCTTGTTGGGCAAAAT
CTCTGCCCTAAAACCTATGCAACCTACTACCTCTTCTAGGAGCCCACTGTATTCTTGCTGGTCTTCTACT
TGGCGGACGCACCTTACATACCTGTTAAGGCAGGCCTTTGAGTTATTTAAAATTAACCACTTAGAATCC
TTTGACTGATATTCTTAAAAGTTCTTGGCCACATGTTCTTGGCAAAAGCCAACAATGTGATCAACTAAAA
GTCAATTTTAAAGTTACTTTAAAATTAATCCAATTACATCAGAGGATTAATAAAGAAACTAACAGTCT
ACCACCTCAGCCTACTGCATACTTTATTTAAAACAATTACCAAACTTCAGATATCAAGTCTCATA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015455](#)

Synonyms: rno-mir-3596b

Locus ID: 100526596



[View online »](#)