

## Product datasheet for **SC403373**

### MIR361 Rat MicroRNA Expression Plasmid (MI0003481)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Rat  
**Vector:** pCMV-MIR (PCMV MIR)  
**Precursor:** GAAGCUUAUCAGAAUCCAGGGGUACUUUUUUUUGAAAAGUCCCCAGGUGUGAUUCUGAUUCGUUUC  
**Mature Sequence:** UUAUCAGAAUCCAGGGGUAC

**Full insert:** >MI0003481 CTGAATGTATGGCAATAAGTCAAGTATCTGAGGGAAAACAATCTTACCATTTGCTTTACATGCCTTGGT  
 TTACAGAGCTTTGCTTAGGCCTATTGCCTCAGTACTTCAGGGCAAGAATGAGGCTAACAGGTGAGTCATG  
 TTACCATCTACCAGACTGCCAGAACTGGGCCTAAATAATAGTCACCTCCCTGTGTTTCAGATGTGCATG  
 TTCTCTCAGTCTGGAGCCCTTTACTTTTCTCACACAACATAGAGATGCAGCAGTAGCAGAGGCTTCACTA  
 TGATTCTTTCTTGGGACTCGGAAGCTTATCAGAATCTCCAGGGGTACTTATTATTTGAAAAGTCCCCAG  
 GTGTGATTCTGATTCTGTTCCCTCTCCTGCTCCTCCTCCTCCTCCTCCTTGGAAATCAGTTTCATCTGA  
 AGCTTTCGCATAAACAGGCACATTTCATTAGCATCTTATCCTGCATAGTTGAGCTTGCATTATTTACTCTC  
 TCGGTTATCTTCATTAATTCTTTCTAATGCTGAAACACCTGAAATCAGTGTCTGCCGGTCACCAAAGGAT  
 CTCTATTAGACTTTGGAAGCAAACAAAGCAGCAAAGACTTCATAAACTTGACTTGAAGGCTAGGTTTATT  
 TCTGGGTTGCTTTGTTTGTGATAATGGGATTGCTTTGA

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)  
**Mammalian Cell Selection:** Neomycin  
**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  
**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.  
**RefSeq:** [MI0003481](#)  
**Synonyms:** rno-mir-361  
**Locus ID:** 100314073



[View online »](#)