

## Product datasheet for **SC403393**

### MIR409a Rat MicroRNA Expression Plasmid (MI0003550)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Rat  
**Vector:** pCMV-MIR (PCMV MIR)  
**Precursor:** GGUACUCGGAGAGAGGUUACCCGAGCAACUUUGCAUCUGGAGGACGAAUGUUGCUCGGUGAACCCUUUUCGGUAUC  
**Mature Sequence:** AGGUUACCCGAGCAACUUUGCAU

**Full insert:** >MI0003550 TCTCCGTTTACAAGGGAGAGTTAGAGCGAGCCGAGTTCTGCACCCACTTCCGCTCCCAGTCAGCTCCTT  
TCTACCTGGGTCTCAAGGTTGACTTCTGGGAAGGCCCTAGATGGTTATTCTGGGAGGTGAAAGGGTGGT  
CGATTGCACTTGGGGAAGTGATGGAATGATTACTCTCACGTGGTGGAGGGTGGCTCCTTGCTGAATGGGG  
AGGTACAAACTCTGGGGGATCTCTGGGATGGTTCTTCTGAAACGCAGCCATCCATCGGCCTTGTCT  
GAGCGTCCAAGCCTCTCCATGGTACTCGGAGAGAGGTTACCCGAGCAACTTTGCATCTGGAGGACGAAATG  
TTGCTCGGTGAACCCCTTTTCGGTATCAAATCCCTCCAGGGAGGCCAATTGGAGGCCAGCACCTCAGGG  
AAGAACGTGGACTCATCTTCGCTCTGGGTATGGGACGGATGGTCGACCAGCTGGAAAGTAATTGTTTCT  
AATGTACTTACCTGGTCCACTAGCCGTGGTGGCCCGCAGCCTGCGCCAGGATGTCACTACCTCAGC  
TCGGCACCTCTGCTGGTACTTGAAGGGAGATCGACCGTGTATATTCGCTTAGCTGACTCGAATAATAC  
ATGTTGATCTTTTCTCAGTATCAATGCTCAGCTTGGAGAAAGCTCT

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)

**Mammalian Cell Selection:** Neomycin

**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.

**RefSeq:** [MI0003550](#)

**Synonyms:** rno-mir-409

**Locus ID:** 100314257



[View online »](#)