

## Product datasheet for **SC403396**

### MIR411 Rat MicroRNA Expression Plasmid (MI0006144)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Rat  
**Vector:** pCMV-MIR (PCMV MIR)  
**Precursor:** CUUGGAGAGAUAGUAGACCGUUAUAGCGUACGCUUUUAUCUGUGACGUAUGUAACACGGUCCACUAACCCUCAGUAUC  
**Mature Sequence:** UAGUAGACCGUUAUAGCGUACG  
**Full insert:** >MI0006144 TGTACCTGAATGCCCTTTCCTTCAGCACACAGAGGGCACGTAAGCAGCATTGTGGAGTAGTAACATCTA  
 CAGCATCGTTTGTCTGCAACCTTTATGAGCTTTTTGGAAGACTTGCCCATTTCCGGGATGAATTGCTGGT  
 GTCAGACCTGACATGAAGATGAGACTCTTGTTTATGTGTGCCTGCCGATGATTGCATGCTCTGTGCC  
 AGCCAGACTCAAGATCCACCAGCTCAATGGCCAGAGAGGGCTTCTGGACGCACCCTCAGCCTTTGTC  
 CTCGGCGTCTCTGTGTGGTACTTGGAGAGATAGTAGACCGTATAGCGTACGCTTTATCTGTGACGTATGT  
 AACACGGTCCACTAACCTCAGTATCAAATCCATCCCCGAGGCCCTGGAGAGGACGCTGTCTTCGGGA  
 GACGCCAGTAGACTCATCATCAGGGCCTCTCCATGACTGGAAGAGATGGTTTTCAATATGGCATGTGT  
 TTTTTGTTTTTTGTTTTTTACTGCATGATGGGAAATGACATAAGCTCCCTCTCCTCTGCTGAGCTC  
 TGCCTTGGGAAGCTTTGCCATCTTTGGGAGGAGGCCATAATGCTTTCTGGGATCAGTGCTCTTCATGAC  
 CTCAGTGCCATGGGTGTCTAGGCAACTTTGGTTAGTAGACTCCAT

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)  
**Mammalian Cell Selection:** Neomycin  
**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  
**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.  
**RefSeq:** [MI0006144](#)  
**Synonyms:** rno-mir-411  
**Locus ID:** 100314107



[View online »](#)