

Product datasheet for **SC403399**

MIR423 Rat MicroRNA Expression Plasmid (MI0006145)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVIR)
Precursor:	GAAGUUAGGCUAGAGGGGCAGAGAGCGAGACUUUUUAUUUUUCCAAAAGCUCGGUCUGAGGCCCCUCAGUCUUGCUUCCU
Mature Sequence:	UGAGGGGCAGAGAGCGAGACUUUU
Full insert:	>MI0006145 CACGTCGGCGTTAGCGTCACGTCGGCGTCGACCACGTCAGCGGAGACGGGAACAAGATGGCGATTCCGG GCAGGCAGTGAGTGATGTGGGAGTTAGGGTCAGGCTGGGGAACAAGAAGCCAGGGAGACCGAGTGTTTAG GCTGAGGCGTCCGAGGTCGCGCGCAGTGACGGGCCATCCCCGCCAGGGTTAGAGGCGGGAGAAGGCAAGG AGTCCCAGCCACATTGATCCCCACCCCGTCTGTATATTTTTCCAGGATGCACGCCCAAGTTTGAG AAACTTGTGAGGAAATAAAGGAAGTTAGGCTGAGGGGCAGAGAGCGAGACTTTTTCTATTTTCCAAAAGCT CGGTCTGAGGCCCTCAGTCTTGCTTCCCTAACCCGCGCTTGAGTTTCTCCCGTTTGGATGCGCCCTGGGG CACCGAGAAGAGAACTCCAGGTTTCTTAAAGGGCCTTGTTCCCGTGGCAGAAGCTCAAGCAAGGCCT GCTTCCCGCCCTCCCTTCTGGTGTAGGAGTGTGGTGCCGAGAAGCCACAGAAAGGGGCAGGCATAAAGG AAACTGCTGTACGCTGCCGTGAAAGACACAAATATTGTTTACCTGCCTGGAAGTTATTCTTGGGCGTGT AAACCTTGCATGCCACGCCCGTGGCCTCTGCATTTTAAAGTCAAAAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0006145
Synonyms:	rno-mir-423
Locus ID:	100314264



[View online »](#)