

Product datasheet for **SC403413**

MIR466b-1 Rat MicroRNA Expression Plasmid (MI0006112)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CAUGUGUAUAUAUGUGUGUGUGUAUGUCCAUGUGUGUAUAUGAAUAUACAUAACACACACACAUACACACACGUGCAAGCACACACA
Mature Sequence:	UAUGUGUGUGUGUAUGUCCAUG
Full insert:	>MI0006112 CACGAATGGCAAACCCGTGGGGAAAGACTCAACAGGATCTGGAGTTCAAACCTCATGCTTCACATCATGG TGAGTTAGAAATCAGACTAGGATTCATGTGATGGCCCTATGTATGTGCATGTGTGTATATGTAAGTGTGT ATGTACATGTGTGTATATGAAGAAACATATACATACACACATACACACACACGCAAGCACACACACACAC ATGAATGGCAAGCACTTGGTAGACAGTATCATGAGGATAATGATTTCAAACCTCAGGCTTCTGGGATTCAT GTGATGTCTCTATATATTTGCATGTGTATATATGTGTGTGTGTATGTCCATGTGTGTATATGAATATACA TACACACACATACACACACGTGCAAGCACACACACAGGAATGGCAAACATTTGGGGGGCAGATTCAGG ATGAGGAGTTAAAGTCATGCTTTACTCATGATTAGTTAGAGCTAGCCTGGGATTCATGTGATTCTCCTGT ATATGTGCATGTGTATATGTCCATGTGTATGTCCATGTCTGTATGTTAATATACATACACAGACATAG AAACACAACCACACATGCTCACACATGCATATGAATGGCAAATACTTGAAGACAGAGCCAACAGTGTCA TGACTTCAAACCTGATTTTTCTGGCATTTCATGGGACACCCCTTTGTATGTGCAT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0006112
Synonyms:	rno-mir-466b-1
Locus ID:	100314093



[View online »](#)