

Product datasheet for **SC403422**

MIR488 Rat MicroRNA Expression Plasmid (MI0006168)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMVIR)
Precursor: AAUCCUCUCUCCAGAUAAUGGCACUCUCAACAAGUUUCUACAUUGUUUGAAAGGCUGUUUCUUGGUCAGAAGACUCU
Mature Sequence: CCCAGAUAAUGGCACUCUCA
Full insert: >MI0006168 C TACTCTGATTCTCTATTCATTTTCAGCTGGTTCTGCCTTCTTGAACCTGAAGGGGACTCCAGCACGGAG
GCAGAAAATGACCCACAGCTGACCTTTTACACCGACCCCTCTCGGAGCCGGAGGCGCAGTCGAGGTAGGG
ATGGTGGGAAGGTGTGTAAAGTGGGGAGCAGGTGGGAAAGGCAGGGTCCACAGCGGTTCTGTGGACTGGA
GCCAECTTCAGCAGGAGGAATGCTTGGATCTCAGCTGGTGCTTTCTCCAGGAAGAATGGAGAGCCTAGG
GTTTCGTTCCGGTAAGAGTGAGAATCCTCTCTCCAGATAATGGCACTCTCAAACAAGTTTCTACATTGTTT
GAAAGGCTGTTTCTTGGTCAGAAGACTCTCAATTTCTTCTGGAACACTCAGCAGGCATGATTTACTGTTC
AGATGAATAAAGTAGATGGGGTTTGATCGATGAATAATATATGTACTCGTCATCAATTCAGCTGCACTGG
CCTTGAGATTGAGTTAGACTGGGATAGGAGGGCATGGTCAGTAAGAATCATCAAGGAGCTGAGGATAACA
TAGTTTTCTTTGGCCTTTAAGTCTTTACTCAGTAGAATCTTGAGACCTAAACCTCCAGACCTTCCCTTC
TCTGCCACAGTTCAGCTAATTTGTTCTCCACACTGACACATCTGTCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0006168](#)
Synonyms: rno-mir-488
Locus ID: 100314120



[View online »](#)