

Product datasheet for **SC403447**

MIR582 Rat MicroRNA Expression Plasmid (MI0012595)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMVMIR)
Precursor: ACUCUUUGGAUACAGUUGUUCAACCAGUUACUAAUCUAAUCUAAUUGUAACCUGUUGAACACUGAACCCAAAGGGUGCAAA
Mature Sequence: UACAGUUGUUCAACCAGUUACU

Full insert:
 >MI0012595 TGTGAAGGTTACCTCCTTCACATTCCTTGAAATATTTTTCTTGTTAAAGTTTCTTTCAGCCAGAGGCCA
 TATATTAATTTCTACTATGTAGGTCAGTCTTTTTCTCTGTTTCTGTCTTTCTCTCCCTCTCTGCCT
 TCCTGTCCGTCATGGTAACCTGCTAATACTAGATGCCAGTGGTTATTTAGCACCTTTCTCTCGGATGTC
 TGCATCTTGTTTCGACGTCATTCTTGACACATTGAAGAGGACAGACTGAGCCCGGCTGAACACCACTGGC
 CAATCTCCCTAGTAATCTGTAATCTTTGGATACAGTTGTTCAACCAGTTACTAATCTAACTAATTGTAAC
 CTGTTGAACAACCTGAACCCAAAGGGTGCAAAGTTGGAAACATTTTCATTGAGAAGTCAGCTGAGGACAAAG
 GAAGCTCGGCCAATCCACTTCAGTGTGGAGAAGGGAAGACTGCTTCTCAAGCACAAACCATGCCACATGG
 TACCTCCGAATGATCGACTGTCCGTTGGTGATGCTGCGTAAGTCTGCATGTATCATCTGTTCTTCCAGTG
 AGGGTTATACATAGTGCAATCCACTAGCCATGAGTGTTTTTTAAACTATCTGAAAATCACAAAAGTTAT
 GAGACTAAAGGATTTGGGAAAATAAGACATATTCCTTTGATTGGTACTTGT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0012595](#)

Synonyms: rno-mir-582

Locus ID: 100314209



[View online »](#)