

Product datasheet for SC403492

MIR711 Rat MicroRNA Expression Plasmid (MI0012621)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMV MIR)
Precursor: CCUAACUUUGAAUCUCUUCGUAGGGUGCUUCAGGAAAGCUGGGACCCUGGGAGAGAUGUAAGUCA
Mature Sequence: GGGACCCUGGGAGAGAUGUAAG
Full insert: >MI0012621 CTTCAAGGCTTCTGCCACTGTGTTTTCTCAGGGTGATCGGGATCCCAGGGCTGGATGGCCGAAATG
GGCTAGATGGGAACCAGGAGCCCTGGCCCTCCAGGGCTACAGGTGAGTGACAGTTTTCTCAGTCCTGC
AGCCCTTCCCTTCACTGGTTAGGTCTTAGTCTTGAGTTCTGGCACCAGCCTGCTCTGTGTCCTTTAG
AGAGCCTTTCTGAGGGGCATATCACTGGTAGGAGGTTTTGGGTTACAGAACTATTGCCTGCCAATGGGC
AGGGACCTCATCCCCTGTGCCTAACTTTGAATCTCTTCGTAGGGTGCTCAGGGAAAGCTGGGGACCT
GGGAGAGATGTAAGTCAAGGGGATGCCAAGGTGGAGGGTGGCTTGGGGATCCAGGAGAGGCCTCATGGG
TGTCACCACTATCAGTGAGAACAGGTCTGCTGGGTCTTTGGAAGGAGAACTCTGTCCCTAGTCAGCAGG
TCACCTGAGTCCCCTTGTGACCTCTGCTTTGCACATCTTTTTTTTTTTTTTTTTTTTTTTTTTTGGTTCT
TTTTTTTGGAGCTGGGGACCGAACCCAGGGCCTTGCTGCTTCCACCACTGAGCTAAATCCCAGCCCTGG
TTTGACATCTAACAGAATGGTTCTGAGGTTCTTCCA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0012621](#)
Synonyms: rno-mir-711
Locus ID: 100314143



[View online »](#)