

## Product datasheet for **SC403494**

### MIR742 Rat MicroRNA Expression Plasmid (MI0006161)

#### Product data:

**Product Type:** MicroRNA Expression Plasmids  
**Species:** Rat  
**Vector:** pCMV-MIR (PCMV MIR)  
**Precursor:** GGCAGUGCUUUACUCACAUGGUUGCUAAUCACGUGAAAUGUAAGUGAAAGCCACCAUGUUGGGUAAAAGUAAAGCUCA  
**Mature Sequence:** UACUCACAUGGUUGCUAAUCA

**Full insert:** >MI0006161 GTGATAACCATATGGAACCATTGTGCTCCATAGAGAAGATAAATGCACTAATTCTGGATAATTCTTCATCTTAATATCTGACGTCATGATTTCAAGTTCTGGAAAGGTTGGGGAACCTTCAAAGTACACTAAGAAACATCTAATAGTGGAACATAAGAGGACAGCATCATTCAGGAATGTGTGGTAAAAGAAAATAATTTATTATTAGAACTTATTTACAAGTTACTTTCAGTGTTCCAGTCTGTTTATTTCAAATGTTGACATGAAGACTGCAGATGTTGTGGTTCTAGAGCGCTGGGCAGTGCTTTACTCACATGGTTGCTAATCACGTGAAATGTAAGTGAAAGCCACCATGTTGGGTAAAGTAAAGCTCATCATGTTCCCTCATTTGTTGTGACCTGGATAGTAAAGTTGAAGAACTACATAGCTTTATTATGACTCCACAGTCTGCTTGCTTGTACTTGAAGAAACCATGGCATCTAAGGGAATTTCTTTTTCTTAATGTTTCAATGTTTCAGAAGTTACCCTTATTCTAGCTTGACAATACCCAAATATACAAATATAAGAACAGAAAATTGAATTCCAATGTTACCAGGATTATGGTTATAAGTGTGAGAAATGAAACCAACCAATGATTTTCATGCCCTACGGAATTAATCAGTTCAACTC

**OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

**E. coli Selection:** Kanamycin (25 ug/mL)  
**Mammalian Cell Selection:** Neomycin  
**Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  
**Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.  
**RefSeq:** [MI0006161](#)  
**Synonyms:** rno-mir-742  
**Locus ID:** 100314116



[View online »](#)