

Product datasheet for **SL500075**

GFP Over-expression Stable Cell Line - HEK293T


Warning: No PDF attribute template selected.

Product data:

_c:	stable-cell-line
_pt:	1
pdf_url:	/datasheet/sl500075.pdf
sds_url:	/sds/sl500075.pdf
Product Name:	GFP Over-expression Stable Cell Line - HEK293T
Description:	GFP Over-expression Stable Cell Line - HEK293T, 1 vial, 3 x 10 ⁶ cells/vial, Cryopreserved
Sizes:	1 vial
Category:	Gene Expression
Product Type:	Stable Cell Lines
selection_marker:	Puromycin (0.8 ug/mL)
Gene Symbol:	GFP
Promoter:	CMV
Parental Cell Line:	HEK293T



[View online »](#)

Product images:

GFP imaging (10x) of HEK293T GFP Stable Cell Lines