

Product datasheet for **SL500077**

GFP Over-expression Stable Cell Line - MDA-MB-231


Warning: No PDF attribute template selected.

Product data:

_c:	stable-cell-line
_pt:	1
pdf_url:	/datasheet/sl500077.pdf
sds_url:	/sds/sl500077.pdf
Product Name:	GFP Over-expression Stable Cell Line - MDA-MB-231
Description:	GFP Over-expression Stable Cell Line - MDA-MB-231, 1 vial, 3 x 10 ⁶ cells/vial, Cryopreserved
Sizes:	1 vial
Category:	Gene Expression
Product Type:	Stable Cell Lines
selection_marker:	Puromycin (0.1 ug/mL)
Gene Symbol:	GFP
Promoter:	CMV
Parental Cell Line:	MDA-MB-231



[View online »](#)

Product images:

GFP imaging (10x) of MDA-MB-231 GFP Stable Cell Lines