

## Product datasheet for **SL500091**

### BCMA and LUC dual-expression stable cell line- HEK293T

#### Product data:

Product Type:	Stable Cell Lines
Gene Symbol:	BCMA
Parental Cell Line:	HEK293T
Promoter:	CMV
Selection Marker:	Hygromycin + Puro (200 ug/ml + 0.8 ug/ml)
Description:	BCMA and LUC dual-expression stable cell line with CAR-T validation - HEK293T, 1 vial, 3 x 10 <sup>6</sup> cells/vial, Cryopreserved



[View online »](#)