


Product datasheet for **SL500093**

Inducible GPR55 over-expression stable cell line- HEK293T

Product data:

Product Type:	Stable Cell Lines
Gene Symbol:	GPR55
Parental Cell Line:	HEK293T
Promoter:	CMV+ Tet-On
Selection Marker:	Hygromycin (200 ug/mL)
Description:	Inducible GPR55 over-expression stable cell line- HEK293T, 1 vial, 3 x 10 ⁶ cells/vial, Cryopreserved

Product images:



Western blotting analysis of induced GPR55-DDK fusion protein in SL500093 stable cells.



[View online »](#)