


## Product datasheet for **SL500095**

### Inducible RHO over-expression stable cell line- HEK293T

#### Product data:

Product Type:	Stable Cell Lines
Gene Symbol:	RHO
Parental Cell Line:	HEK293T
Promoter:	CMV+ Tet-On
Selection Marker:	Hygromycin (200 ug/mL)
Description:	Inducible RHO over-expression stable cell line- HEK293T, 1 vial, 3 x 10 <sup>6</sup> cells/vial, Cryopreserved

#### Product images:



Western blotting analysis of induced RHO-DDK fusion protein in SL500095 stable cells with an anti-DDK antibody.



[View online »](#)