


Product datasheet for **SL500096**

Inducible OTUB1 over-expression stable cell line- HEK293T

Product data:

Product Type:	Stable Cell Lines
Gene Symbol:	OTUB1
Parental Cell Line:	HEK293T
Promoter:	CMV+ Tet-On
Selection Marker:	Hygromycin (200 ug/mL)
Description:	Inducible OTUB1 over-expression stable cell line- HEK293T, 1 vial, 3 x 10 ⁶ cells/vial, Cryopreserved

Product images:



Western blotting analysis of induced OTUB1-DDK fusion protein in SL500096 stable cells.



[View online »](#)