

## Datasheet for Neuro-2a/Cas9-Rosa26-Puro Safe Harbor Cell Line

**Catalog number:** SL550006

**Product:** Mouse Neuro-2A cell line stably expressing CRISPR Cas9 nuclease from ROSA26 locus.


**Description:** This product is a cell line stably expressing the CRISPR Cas9 nuclease. Cas9 is integrated at the mouse Rosa26 Safe Harbor locus. This cell line also expresses the puromycin resistance gene. In combination with separately transfected or transduced single guide RNAs (sgRNAs), this cell line will sustain double-strand DNA breaks (DSBs) at targeted genome sites. Though we could not detect Cas9 expression by Western blot, we found the basal Cas9 activity by T7 Endonuclease I assay when sgRNA expression at high levels using plasmid transfection. This cell line can be used *in vitro* for gene knockout, transgene knockin, mutagenesis, transgene integration, or other genome editing-related applications.

**Quantity:** 1 vial of  $2 \times 10^6$  cells; frozen

**Shipping conditions:** Dry ice

**Storage conditions:** Liquid nitrogen vapor phase. Remove the item from the dry ice packaging and check all items for damage and leakage. Place immediately into storage at or below  $-140^\circ\text{C}$ , preferably into the liquid nitrogen vapor phase, until use.

**Transgene integration:**



**Source of parental line:**

Neuro-2a  
 Organism: *Mus musculus*, mouse  
 Tissue: B r a i n  
 Cell type: Neuroblast

- Quality control:** >95% viability before freezing. All cells were tested and found to be free of mycoplasma, bacteria, viruses, and other toxins.
- Safety instructions:** To ensure safety, protective gloves, clothing, and a face mask should be worn when handling frozen vials. Some leakage may occur into the vial during storage. The liquid nitrogen will be converted to gas upon thawing. Due to the nature of nitrogen gas, pressure may build within the vial and possibly result in the vial exploding or losing its cap. This may cause flying debris.
- Thawing procedure:** The vial of cells should be thawed in a 37 °C water bath with gentle agitation. For optimal performance, the vial should be thawed in under two minutes. Ensure that the cap of the vial did not loosen upon thawing, and re-tighten as needed. Spray the vial with 70% EtOH and wipe off. Repeat. Using aseptic technique, add the contents of the vial to 9 ml of complete growth medium (without selection). Centrifuge for 5 min. at 125 x g. Aspirate the medium, being careful not to disturb the pellet. Resuspend in 10 mL of complete growth medium, and place into a culture vessel of your choice. Only add selection to the medium after 24 hours in culture.

**Culture conditions:**

**Complete Growth Medium**

The base medium for this cell line is DMEM. For optimal growth and maintenance of selection, add the following components to the base medium: fetal bovine serum to a final concentration of 10%.

**Selection**

Puromycin to a final concentration of 1µg/mL.

**Culture temperature:**

37 °C with 5% CO<sub>2</sub>

**Subculture:**

Rinse the cells with PBS without cations, digest cells with 0.25% (w/v) Trypsin-EDTA (0.53 mM) solution and split at 1:3 to 1:10 ratio.

**Cryopreservation:** Freeze slowly in complete growth medium supplemented with 5% (v/v) DMSO.

**Mycoplasma:** Negative  
(MycoAlert Mycoplasma Detection Kit from Lonza)

---

## Product QC:

### 1. Junctional PCR (to confirm the Cas9 gene integration into AAVS1 site)


#### (1) 5' Junctional PCR

One primer recognizes the region of the chromosome outside of the 5' homology arm region, while the other primer recognizes the Cas9-plasmid region to confirm the ROSA26 site integration of the Cas9 expressing cassette

Junction-PCR 5'F: AGGGAGCGGAAAAGTCTCCA

Junction-PCR 5'R: GCGTACTTGGCATATGAT

**Predict product length: 1129 bp**




#### (2) 3' Junctional PCR to confirm 3' integration

one primer from chromosomal outside of the 3' homology arm region, the other primer from the Cas9-plasmid region to confirm the AAVS1 site integration of Cas9 expressing cassette

Junction-PCR 3'F: GGTGCCTGAAATCAACCTCTGGA

Junction-PCR 3'R: TCAAGCCAGTCCAAGAGAAAGCA


**Predict product length: 1604 bp**



## 2. T7 Endonuclease I (T7 E1) Assay

### (1) P2X2r Site T7 E1 Assay

sgRNA targeting to P2X2r gene was transduced into Neuro-2a/Cas9-Rosa26-Puro Safe Harbor cell line by transient transfection. After transfection, Cas9 protein was induced for 72 hours. P2X2r gene was cut by CAS9 expressed inside the cells and repaired through NHEJ with mutation. A 730 bp P2X2r gene fragment from PCR was then tested by T7 Endonuclease I (T7 E1) Assay. The T7 E1 cleavage will results in two additional bands: one ~494 bp and the other ~236 bp.



**Citation of product:** If use of this item results in a publication, please use this information:  
CRISPR Cas9 stable Neuro-2a/Cas9-Rosa26-Puro Safe Harbor cell line (SL550006,  
OriGene Technologies, Inc., Rockville, MD).

#### Limited Use License

A limited use license is granted to the Buyer of the Product. The Product shall be used by the Buyer for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use. The Product must not be resold, repackaged or modified for resale, or used to manufacture commercial products without prior written consent from OriGene. This Product should be used in accordance with NIH guidelines developed for recombinant DNA and genetic research. Use of any part of the Product constitutes acceptance of the above terms.