

Product datasheet for **SM1001P**

Human IgA (Secretory component) Mouse Monoclonal Antibody [Clone ID: MC29/12]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (Secretory component) Mouse Monoclonal Antibody [Clone ID: MC29/12]
Clone Name:	MC29/12
Applications:	ELISA, IHC
Recommended Dilution:	Immunohistochemistry on frozen and paraffin sections. ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	Spleen cells from immunised BALB/c mice were fused with cells of the NS-1 myeloma cell line.
Isotype:	IgG1
Formulation:	Tris buffered saline containing 0.09% sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)