

## Product datasheet for **SM1008P**

### Human IgG1 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 8c/6-39 (HP6091)]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG1 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 8c/6-39 (HP6091)]
Clone Name:	8c/6-39 (HP6091)
Applications:	ELISA, ID, IHC
Recommended Dilution:	<b>ELISA.</b> <b>Immunodiffusion.</b> <b>Immunohistochemistry on Frozen Sections.</b> <b>Not suitable for Paraffin Sections.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified Human IgG1
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)