

## Product datasheet for **SM1010P**

### Human IgG2 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: HP6014]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG2 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: HP6014]
Clone Name:	HP6014
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry on Frozen Sections.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified Human IgG2.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)