

Product datasheet for **SM1017F**

Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-2]
Clone Name:	LO-MG1-2
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 10 µl of 5 µg/ml diluted antibody to label 10 ⁶ cells in 100 µl.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified Murine IgG1 from BALB/c mice
Isotype:	IgG1
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction from Tissue Culture Supernatant Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1 Absorption emission: 490 nm / 525 nm
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use. Avoid repeated freezing and thawing.



[View online »](#)