

Product datasheet for **SM1017PX**

Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-2]
Clone Name:	LO-MG1-2
Applications:	ELISA, IHC
Recommended Dilution:	Useful as second antibody in immunoassays, such as ELISA's (good binding on plastic), Immunohistochemistry on Frozen Sections and for Affinity Chromatography of murine IgG1.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)