

Product datasheet for **SM1018HRP**

Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-7]
Clone Name:	LO-MG2a-7
Applications:	ELISA
Recommended Dilution:	ELISA: 0.5 µg/ml as coating antibody.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG from BALB/c mice. Spleen cells from immunised LOU/c rats were fused with cells of the rat IR983F myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.01% Thiomersal and 50% Glycerol Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)