

## Product datasheet for **SM1018P**

### Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-7]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-7]
Clone Name:	LO-MG2a-7
Applications:	ELISA
Recommended Dilution:	ELISA: 5µg/ml as coating antibody.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG from BALB/c mice.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography from tissue culture supernatant.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)