

Product datasheet for **SM1021B**

Mouse IgG2b (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2b-2]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Mouse IgG2b (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2b-2] |
| Clone Name: | LO-MG2b-2 |
| Applications: | ELISA |
| Recommended Dilution: | ELISA: 500 ng/ml as a detection antibody. |
| Reactivity: | Mouse |
| Host: | Rat |
| Immunogen: | Purified IgG from BALB/c mice. Spleen cells from immunised LOU/c rats were fused with cells of the rat IR983F myeloma cell line. |
| Isotype: | IgG1 |
| Formulation: | PBS containing 0.09% Sodium Azide and 50% Glycerol Label: Biotin State: Liquid purified IgG |
| Concentration: | lot specific |
| Purification: | Affinity chromatography on Protein A |
| Conjugation: | Biotin |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)