

Product datasheet for **SM1021P**

Mouse IgG2b (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2b-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2b-2]
Clone Name:	LO-MG2b-2
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml as a detection antibody.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG from BALB/c mice. Spleen cells from immunised LOU/c rats were fused with cells of the rat IR983F myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)