

Product datasheet for **SM1022B**

Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-7]
Clone Name:	LO-MG3-7
Applications:	ELISA
Recommended Dilution:	ELISA: 0.5-1.0 µg/ml.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Polyclonal and Monoclonal Mouse IgG3.
Isotype:	IgM
Formulation:	PBS Label: Biotin State: Liquid purified IgM fraction from Tissue Culture Supernatant Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)