

## Product datasheet for **SM1022F**

### Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-7]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-7]
Clone Name:	LO-MG3-7
Applications:	ELISA, IHC
Recommended Dilution:	Useful as second antibody in <b>Immunoassays</b> , such as <b>ELISA's</b> (good binding on plastic) and <b>Immunohistochemistry on Frozen Sections</b> . The antibody can be used at 5 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgM
Formulation:	PBS Label: FITC State: Purified State: Liquid purified Ig fraction Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide Molar ratio: 4.6
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)