

Product datasheet for **SM1023B**

Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]
Clone Name:	LO-ME-3
Applications:	ELISA
Recommended Dilution:	ELISA (Use at 0.5 µg/ml as a detection antibody, avidity 3 x 10e9 M-1).
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified mouse IgE
Isotype:	IgG1
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction from Tissue Culture Supernatant Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)