

Product datasheet for **SM1023F**

Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]
Clone Name:	LO-ME-3
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: Use 50µl of 5µg/ml diluted antibody to label 10e6 cells in 100µl.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified Mouse IgE from BALB/c mice
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide and 50% Glycerol Label: FITC State: Liquid purified IgG fraction from Tissue Culture Supernatant Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)