

Product datasheet for **SM1023FX**

Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]
Clone Name:	LO-ME-3
Applications:	ELISA
Recommended Dilution:	Useful as second antibody in immunoassays , such as ELISA's
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS Label: FITC State: Liquid purified Ig fraction Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C. This products is photosensitive and should be protected from light.



[View online »](#)