

Product datasheet for **SM1023HRP**

Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]
Clone Name:	LO-ME-3
Applications:	ELISA
Recommended Dilution:	ELISA: 500ng/ml.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified Mouse IgE. Spleen cells from immunised LOU/c rats were fused with cells of the rat IR983F myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 50% Glycerol and 0.01% Thiomersal. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)