

Product datasheet for **SM1023PX**

Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-ME-3]
Clone Name:	LO-ME-3
Applications:	ELISA, IHC
Recommended Dilution:	Useful as second antibody in Immunoassays , such as ELISA's (good binding on plastic) and Immunohistochemistry on Frozen Sections .
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	State: Purified State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)