

Product datasheet for **SM1025P**

Mouse IgM (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MM-9]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MM-9]
Clone Name:	LO-MM-9
Applications:	ELISA
Recommended Dilution:	ELISA: 5ug/ml (as a coating antibody). Immunohistochemistry on frozen sections.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Polyclonal and monoclonal mouse IgM. Spleen cells from immunised LOU/c rats were fused with cells of the IR983F rat myeloma cell line.
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on mouse IgM/Sepharose 4B
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)