

Product datasheet for **SM1026AP**

Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]
Clone Name:	HP6050
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/500-1/1,000. Western blotting. Dot- and Slot-Immunoblotting Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 Label: AP State: Liquid IgG fraction Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)