

Product datasheet for **SM1026B**

Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]
Clone Name:	HP6050
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	ELISA: 1/5,000-1/10,000. Western blotting. Dot- and Slot-Immunoblotting. Immunocytochemistry. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)